

Quartz tuning
fork vibration
amplitude as a
limitation of
spatial resolution
of shear force
microscopes

J.-M Friedt, É.
Carry, Z. Sadani,
B. Serio, M.
Wilm, S.
Ballandras
FEMTO-ST

Quartz tuning fork vibration amplitude as a limitation of spatial resolution of shear force microscopes

Finite Element
Analysis

Experimental
results

Data processing

Stroboscopic
method

Conclusion

J.-M Friedt, É. Carry, Z. Sadani, B. Serio, M. Wilm, S. Ballandras
FEMTO-ST

23 mars 2005

Why ?

Quartz tuning
fork vibration
amplitude as a
limitation of
spatial resolution
of shear force
microscopes

J.-M Friedt, É.
Carry, Z. Sadani,
B. Serio, M.
Wilm, S.
Ballandras
FEMTO-ST

Finite Element
Analysis

Experimental
results

Data processing

Stroboscopic
method

Conclusion

- SPM usually use the physical quantity under investigation as probe-distance indication
- this is fine on homogeneous surface (constant physical quantity)
- shear-force microscopy uses a resonator for independent probe-distance feedback
→ usable for a wide range of applications (SNOM, SECM, STM ...)

Quartz tuning
fork vibration
amplitude as a
limitation of
spatial resolution
of shear force
microscopes

J.-M Friedt, É.
Carry, Z. Sadani,
B. Serio, M.
Wilm, S.
Ballandras
FEMTO-ST

Finite Element
Analysis

Experimental
results

Data processing



Stroboscopic
method

Conclusion

Shear force microscopy

Feedback on one property of a quartz resonator (current magnitude or phase) to keep the probe-surface distance constant : the resonator is disturbed by the forces acting on the tip

⇒ modification of the transfer function of the resonator. The feedback signal (probe-distance) is recorded for topography monitoring.



But ...

Shear force microscopy has not displayed the excellent resolution of other scanning probe microscopies

⇒ requires a good understanding of the behavior of the probe and its interaction with the surface

Quartz tuning fork vibration amplitude as a limitation of spatial resolution of shear force microscopes

J.-M Friedt, É. Carry, Z. Sadani, B. Serio, M. Wilm, S. Ballandras
FEMTO-ST



Finite Element Analysis

Experimental results

Data processing

Stroboscopic method

Conclusion



- size of the probe ?
- “leakage” of the near field (evanescent) physical property ?
- vibration amplitude of the probe ?

K. H. Choi, J.-M Friedt, F. Frederix, ... *Simultaneous Atomic Force Microscope and Quartz Crystal Microbalance Measurement Applied*

Physics Letters (Vol 81, No 7, 12 Aug 2002)

Quartz tuning
fork vibration
amplitude as a
limitation of
spatial resolution
of shear force
microscopes

J.-M Friedt, É.
Carry, Z. Sadani,
B. Serio, M.
Wilm, S.
Ballandras
FEMTO-ST

Finite Element
Analysis


Experimental
results

Data processing

Stroboscopic
method

Conclusion

Modulef based dynamic simulation : free tuning fork




Quartz tuning
fork vibration
amplitude as a
limitation of
spatial resolution
of shear force
microscopes

J.-M Friedt, É.
Carry, Z. Sadani,
B. Serio, M.
Wilm, S.
Ballandras
FEMTO-ST

Finite Element
Analysis
Experimental
results
Data processing
Stroboscopic
method
Conclusion

Modulef based dynamic simulation : loaded free tuning fork




Quartz tuning fork vibration amplitude as a limitation of spatial resolution of shear force microscopes

J.-M Friedt, É. Carry, Z. Sadani, B. Serio, M. Wilm, S. Ballandras FEMTO-ST

Finite Element Analysis
 Experimental results
 Data processing
 Stroboscopic method
 Conclusion

Interferometric methods



Quartz tuning
fork vibration
amplitude as a
limitation of
spatial resolution
of shear force
microscopes

J.-M Friedt, É.
Carry, Z. Sadani,
B. Serio, M.
Wilm, S.
Ballandras
FEMTO-ST


Finite Element
Analysis

Experimental
results


Data processing

Stroboscopic
method

Conclusion



Raw data





Quartz tuning
fork vibration
amplitude as a
limitation of
spatial resolution
of shear force
microscopes

J.-M Friedt, É.
Carry, Z. Sadani,
B. Serio, M.
Wilm, S.
Ballandras
FEMTO-ST

Finite Element
Analysis
Experimental
results
Data processing
Stroboscopic
method
Conclusion

Linking model and experimental data



abscissa is graduated in voltage from 400 to 9000 mV (experimental data), which is also equal to (simulated data) a vibration amplitude of $\lambda/21 = 23$ nm to $\lambda/1.7 = 290$ nm ($\lambda = 488$ nm in this experiment).

Quartz tuning fork vibration amplitude as a limitation of spatial resolution of shear force microscopes

J.-M Friedt, É. Carry, Z. Sadani, B. Serio, M. Wilm, S. Ballandras FEMTO-ST

Finite Element Analysis


Experimental results

Data processing

Stroboscopic method

Conclusion

Loaded tuning fork







abscissa spans from 100 mV to 7600 mV amplitude (experiment) which is also equal to $\lambda/126=4$ nm to $\lambda/5=97$ nm (here $\lambda = 488$ nm).

Quartz tuning fork vibration amplitude as a limitation of spatial resolution of shear force microscopes

J.-M Friedt, É. Carry, Z. Sadani, B. Serio, M. Wilm, S. Ballandras FEMTO-ST

Stroboscopic method

- 1 acquire images of the moving surface phase-synchronized with the driving voltage
- 2 oversample each line and intercorrelate images
- 3 look for the maximum of intercorrelation and find the best sine-wave fit
- 4 repeat for each line of the image



Results : $0.5 V_{pp}$ – $1.8 V_{pp}$ – $5 V_{pp}$ → displacement amplitude
 350 nm–850 nm–3000 nm ($Q = 4500$).

Conclusion and perspectives

- we have developed the basic Finite Element Model of a tip-loaded tuning fork
- we have experimentally measured the vibration amplitude of a tuning fork

Further developments include :

- measuring the vibration amplitude as a function of probe-surface distance
- adding an external force acting on the tip of the probe to our model
- experimentally observe possible spatial resolution loss at high driving-voltage amplitude
- is the tuning fork usable as a scanner ? (for an $N \times N$ pixel image, we must sample at $N \times 32768$ Hz to get a framerate of $32768/N$ image/s : $N = 128 \Rightarrow 8.5$ Msamples/s and 512 fps!)

